



**REF K100-0104**

## **E-Z Chek<sup>®</sup> 0 - 14 pH Test Strips**

### **Instructions for Use**

#### **Indications for Use**

K100-0104 E-Z Chek<sup>®</sup> 0 - 14 pH Test Strips are indicated for use in testing the pH of water used for preparing dialysate. There are no other indications for use.

#### **Directions**

- 1) Remove one strip from the bottle and immediately recap the bottle.
- 2) Dip the strip into the test solution until the color on the indicator pads no longer changes.
- 3) Remove the strip and compare it to the color chart on the bottle.
- 4) Match the strip as closely as possible to one of the color scales on the bottle. Read the pH listed above the matched colors. Colors that are between two scales indicate a pH that is between the two values for those scales.

#### **Warnings and Cautions**


- Do not touch the indicator pads.
- Keep all unused test strips in the original bottle.
- Replace cap immediately/tightly after removing a strip.
- Do not use test strips from an opened or unopened bottle after expiration date printed on the bottle label.
- Do not allow the test strip to come in contact with liquids or with work surfaces that may contaminate the strip.
- Do not leave test strips in areas exposed to vapors of any type.
- *Caution:* When used as a medical device, Federal Law restricts this device to sale by or on order of a physician.

#### **Interferences**

No interferences are known.

#### **Storage**

- E-Z Chek<sup>®</sup> 0 – 14pH test strips must be kept in the original bottle with the lid tightly closed.

39°F (4°C)  86°F (30°C)



- Do not use a test strip from an opened or unopened bottle after the expiration date.
- Lot number and expiration date are printed on the bottle label.

**RPC • Tucson, AZ • Minneapolis, MN • USA • (520) 888-5551 • [orders@rpc-rabrenco.com](mailto:orders@rpc-rabrenco.com)**

**Product Quality at Discount Pricing<sup>SM</sup>**