

Micro-X® Potency & Residual Test Strips

RPC's Micro-X[®] test strips are designed for use in testing Micro-X[®] cold sterilant and other substantially equivalent peracetic acid (PAA) germicides such as Renalin, Minncare, Peracidin, and Puristeril. Use the K100-0105 potency test strips to make sure the proper concentration of the PAA germicide is used when disinfecting dialyzers, water systems, dialysis machines, etc. After the disinfection process and minimum exposure times are complete, the K100-0100 Micro-X[®] residual test strips are used to test the rinse solution to determine if the residual levels of PAA germicide have been rinsed to a safe level (per the germicide Instructions for Use).

- Compatible with all dialysis peracetic acid germicides
- Five color blocks for easy, accurate testing of residual PAA germicide
- Quick and easy "dip and read" tests with results in 10 seconds or less



- Save time and money on rinsing as residual strips "turn off" below 1 ppm.
- Accurate potency strips turn positive at or above proper PAA strength
- FDA cleared for testing of PAA germicides used in dialysis.

Micro-X[®] Potency & Residual Test Strips are supplied in easy to use containers designed for use with RPC's popular CapKeeper[®] holder. A Certi-Chek[™] certificate of conformance for QC validation is available for each test strip type, by lot number. Labeled product shelf life is valid for unopened containers and for opened containers that are immediately resealed and kept in storage conditions as stated in the instructions for use. RPC test strips provide you with multiple options for measuring and monitoring critical dialysis fluids. RPC test strips help maximize quality in your dialysis program and minimize risk when testing dialysis fluids.

For Orders, Technical Support, or Product Catalog Call: 1-800-647-3873, Fax: 877-352-5557 or Email: orders@rpc-rabrenco.com

Savings, solutions, and satisfaction [™] all from RPC, the leader in dialysis test strip applications.		
Part Number	Description	Packaging
K100-0100	Micro-X [®] Residual Test Strips	100 per bottle
K100-0105	Micro-X [®] Potency Test Strips	100 per bottle