




CapKeeper® Test Strip Cap & Container Organizer Instructions For Use



Push/pull-off Cap Style

1. Push capped test strip container completely into Capkeeper® cap cover.
2. To open/close strip container: grip and squeeze cap cover while pulling/ pushing on test strip container.

Rotate-on/off Cap Style

1. Push capped test strip container completely into CapKeeper® cap cover.
2. To open/close strip container: grip and squeeze cap cover while rotating test strip container.

RPC • Tucson, AZ • Minneapolis, MN • USA • (520) 888-5551 • orders@rpcrabrenco.com

Product Quality at Discount PricingSM