



## Citra-Quik™ Dry Powder Citric Acid For Cleaning Dialysis Machine Fluid Pathways

Bicarbonate-based dialysate results in carbonate scaling within hemodialysis machines over time. This precipitate will eventually cause the dialysis machine to malfunction, resulting in problems associated with downtime. Historically, maintenance personnel have used various chemical solutions to remove scale deposits. The following list of key benefits illustrates why only citric acid from RPC should be used.

**SAVE!**

- ◆ Replaces hazardous chemicals such as nitric and oxalic acid.
- ◆ No obnoxious fumes as with vinegar.
- ◆ Dissolves quickly, easily in AAMI Quality Water.



- ◆ Faster and better at removing deposits and won't harm equipment.
- ◆ Costs 20-50% less than vinegar or other cleaners.
- ◆ Save on storage space and shipping costs.

**IMPROVE!**

Citric acid from RPC is also flexible. It allows mixing to the most effective concentration (typically 5% - 50%). Concentration with commercial pre-mixed liquid solutions is fixed. An informal industry leading study recently showed that RPC's citric acid went totally into solution many times faster than competitor's citric acid, even at cold water temperatures. This can save a lot of time and prevent sucking up citric acid particulate into the dialysis machines. Citric acid is excellent for other purposes also, such as removing corrosion from stainless steel parts or dried dialysate from plastic parts.

Flat packaging and dry powder save on storage space and shipping costs compared to pre-mixed solutions. Each bag contains the highest quality USP grade citric acid which is made in the USA via controlled quality processing. It is packaged to keep out moisture and prevent leakage.

**For Product Catalog, Orders or Technical Support:  
Call: 1-800-647-3873, or Fax: 1-877-352-5557**

*Available in five convenient package sizes:*

P/N J100-0190	Citra-Quik™ Powder 190 gram bag, 30 per case
P/N J100-0408	Citra-Quik™ Powder 408 gram bag, 30 per case
P/N J100-0747	Citra-Quik™ Powder 747 gram bag, 20 per case
P/N J100 1000	Citra-Quik™ Powder 1000 gram (1 kg) bag, 15 per case
P/N J100 1050	Citra-Quik™ Powder 50 lb. bag (boxed)

*For Free Technical Support Information visit RPC's Web site:  
[www.rpc-rabrenco.com](http://www.rpc-rabrenco.com)*

© 2017 By RPC