

CathGuard[™] Patient Catheter Protective Pocket

CathGuard[™] is a unique, water-resistant, soft fabric, protective pocket. It provides maximum patient comfort while protecting and storing the exterior lumen of the patient's central line catheter between dialysis treatments. A patient's catheter can be quickly and easily placed into the pocket. Each CathGuard[™] pocket has its own self-adhesive strips for attaching the pocket to the patient's skin. CathGuard[™] adhesive strips are designed to attach the pocket securely to the patient, yet make the pocket easy to remove without causing discomfort for the patient.

- Protects patient and catheter. Soft, waterresistant fabric pocket.
- Eliminates mess and fuss of dealing with tape and gauze.
- Adhesive strips secure pocket to patient without discomfort on removal.



- Saves staff time by simplifying patient offdialysis procedure.
- Easy-to-use: insert catheter; remove backing; attach to patient's skin.
- Cost-effective, unique product. Well-liked by patients and staff.

Use of CathGuard[™] helps to immobilize the patient's catheter and may aid in reducing trauma to the catheter surgical insertion site. In addition, CathGuard[™] keeps the catheter clean, may help prevent infections, safeguards the integrity of the catheter caps, and presents a flatter, less-visible profile.

RPC also offers high quality, cost-effective non-injectable catheter caps (P/N 7001-3915, 100 per case).

For Product Catalog, Orders or Technical Support:

Call: 1-800-647-3873, Fax:1-877-352-5557 or Email: orders@rpc-rabrenco.com

