




E-Z Grip® Port Caps & Renatron Reprocessing Connector

Dialysis reuse managers asked us at RPC to design and manufacture dialyzer port caps that were easier to use and reuse. We listened and are pleased to offer our E-Z Grip® Port Caps. A partial list of the many E-Z Grip® benefits is shown below. In addition, RPC offers an **O-ring free** reprocessing connector for use with the Renatron. This connector can be used with or without the RPC Sani-lok® Barrier Adapter.

E-Z GRIP® DIALYZER PORT CAPS

- ◆ Easy to grip, rigid plastic for fast installation & removal.
- ◆ Compatible with heat and germicide reuse.
- ◆ No phthalate esters such as DEHP plasticizer.
- ◆ Caps & connector have no O-rings to leak, or fall out of place.
- ◆ Save money compared to Fresenius caps with O-rings.

Ask for Free Samples!



<u>Part Number</u>	<u>Description</u>	<u>Packaging</u>
7000-1051	Reprocessing Connector (no o-ring, replaces Renatron P/N 78399-699)	50 per bag
7000-0009	E-Z Grip® Dialysate Port Cap Twist Lk-Rigid (Fresenius dialyzers only)	500 per bag

Zebra Printer Labels, Ribbon & Cleaning Card

Renalog Labels

Used in conjunction with Item# 7000-2403 and Zebra TLP Label Printer. Labels are easy to peel off roll and will adhere securely to all dialyzers.



Ribbon Printer

Used in conjunction with Item# 7000-2102 and Zebra TLP 2844 Label Printer these black printer ribbons will display patient information and reprocessing cycles with consistent quality and legibility.



- ◆ Save money; low cost, high quality
- ◆ Equal to or better than manufacturer's versions

<u>Part Number</u>	<u>Description</u>	<u>Packaging</u>
7000-2102	Renalog Labels for Zebra TLP 2844 label printer - 1750 labels per roll	1 roll
7000-2403	Black Printer Ribbon for Zebra TLP 2844	1 ribbon