



For Free Technical Support
Information Visit RPC's Web-site:
www.rpc-rabrenco.com

E-Z Chek® Sensitive Total Chlorine & Chloramines Test Strips

E-Z Chek® Sensitive Total Chlorine & Chloramines Test Strips provide a quick and convenient means for indicating low levels of total chlorine in water. Monitoring between carbon tanks for chlorine/chloramines break through is the primary application for the E-Z Chek® Sensitive. This strip is a low cost option for go/no go testing of carbon bed performance (color block at 0.1 ppm limit). **A second application** for these strips is to provide a convenient means for indicating the concentration of residual chlorine (0.5 ppm chlorine limit) following bleach disinfection of hemodialysis equipment and water distribution loops. **A third application** for the E-Z Chek® Sensitive strips is for use in monitoring chlorine levels in the raw water delivered to your facility's water room. Chlorine is used by municipalities to disinfect the raw water supply. The EPA raw water maximum allowable concentration is 4 ppm chlorine. This level should be monitored periodically.

- ◆ **Eliminate messy and time consuming reagent tests**
- ◆ **Four color blocks for easy, accurate testing of total chlorine**
- ◆ **Verify compliance with recommendations for residual chlorine in rinse and raw water.**



- ◆ **Easy “dip and read” test procedure. Results in only 40 seconds**
- ◆ **Accurate, low cost strips detect levels at less than 0.1 ppm**
- ◆ **Better color contrast between test values when compared with competitors strips**

The E-Z Chek® Sensitive Total Chlorine & Chloramines Test Strips are supplied in a ready to use form. A certificate of conformance for QC validation is available for each test strip type, by lot number. When using RPC test strips all test strip instructions for use, including storage and handling requirements, must be followed. RPC test strips provide you with multiple options to measure and monitor the critical testing, needed to ensure the highest quality dialysis program.

For Orders, Technical Support, or Product Catalog

Call: 1-800-647-3873, Fax: 877-352-5557 or Email: orders@rpc-rabrenco.com

Savings, solutions, and satisfactionSM ... all from RPC, the leader in dialysis test strip applications.

Part Number	Description	Packaging
K100-0106	Total Chlorine Test Strips, E-Z Chek® Sensitive	100 per btl.

Product Quality at Discount PricingSM