



For Free Technical Support
Information Visit RPC's Web-site:
www.rpc-rabrenco.com

E-Z Chek® 6.8–8.5 pH Dialysate & Concentrate Test Strips

RPC's unique *E-Z Chek® 6.8 - 8.5 pH Test Strips* are an accurate, easy, and cost-effective means for validating the pH of dialysate against the AAMI (RD52) specified range of 6.9 to 7.6. *Each test strip has two color pads* designed to give maximum color differentiation and clarity between the six pH increments shown on the bottle color scale for the K100-0117 and four pH increments on the K100-0117CT. On the color scale, a red line square surrounds the 6.8 and 7.8 pH values to indicate "no go" values outside the AAMI recommended range for dialysate. The color scale also includes two bicarbonate concentrate color blocks at 8.0 and 8.5 pH. RPC's E-Z Chek® pH Test Strips are packaged 100 strips in a convenient, round, flip-top plastic bottle rather than an awkward-to-use plastic flat pack. For optimal accuracy, two versions are available. One version is for use with bicarb/acetate dialysate (K100-0117) and the other is for use with bicarb/citrate dialysate (K100-0117CT).

- ◆ **Fast 10 sec. "dip and read" test for pH in dialysate.**
- ◆ **Color values designed for easy, accurate testing of 6.8 - 7.8 pH levels.**
- ◆ **Verify compliance with AAMI recommendation for dialysate pH.**



- ◆ **Much more accurate than competing pH test strips.**
- ◆ **FDA cleared for testing dialysate/concentrate pH.**
- ◆ **Bicarb/acetate and bicarb/citrate versions available for optimal accuracy and precision.**

E-Z Chek® 6.8-8.5 pH Dialysate & Concentrate Test Strips are supplied in a ready to use form. A Certi-Chek™ certificate of conformance for QC validation is available for each test strip type, by lot number. When using RPC test strips all test strip instructions for use, including storage and handling requirements, must be followed. RPC test strips give you multiple options for measuring and monitoring critical dialysis fluids. RPC test strips help maximize quality in your dialysis program and minimize risk when testing dialysis fluids.

For Orders, Technical Support, or Product Catalog

Call: 800-647-3873, Fax: 877-352-5557 or Email: orders@rpc-rabrenco.com

Part Number	Description	Packaging
K100-0117	E-Z Chek® 6.8 – 8.5 pH Test Strips (bicarb/acetate dialysate)	100 per bottle
K100-0117CT	E-Z Chek® 6.8 – 8.5 pH Test Strips (bicarb/citrate dialysate)	100 per bottle

Product Quality at Discount PricingSM