



For Free Technical Support Information, Visit RPC's Web-site: [www.rpc-rabrenco.com](http://www.rpc-rabrenco.com)

## Micro-X® Dialyzer Reprocessing Cold Sterilant

Simplify life, save time and save money, with one-stop shopping and unequalled technical support for your dialyzer reprocessing program. By adding Micro-X® to its product line, Reprocessing Products Corporation (RPC) now offers you far more reprocessing resources than any other company. Dialysis professionals recognize RPC as the number one company for giving customers exceptional service and a high quality family of cost-effective products for dialyzer reprocessing. Rely on RPC: use Micro-X® in your reprocessing program.

◆ **One-stop shopping for reprocessing, water, & other technical supplies.**

◆ **Same proven peracetic acid formula in use for many years.**

◆ **Expert support from the “gurus”... all methods of dialyzer reprocessing.**



◆ **Fast (3 - 6 day) shipments to all areas of the United States.**

◆ **FDA clearance for automated & manual reprocessing systems.**

◆ **No need to change procedures or protocols currently used.**

Minimize procedural changes in your dialysis center: Micro-X® has identical chemical properties and intended use as Renalin®. It also has US Food & Drug Administration marketing clearance for use in all manual or automated reprocessing systems, e.g., Renatron, Renatron II, Echo, DRS-4, etc.

RPC offers a family of the highest quality reprocessing products, including port caps designed specially for reprocessed dialyzers; Auto-Vent® Port Cap; Sani-lok® Barrier Adapter; Micro-X® Potency & Residual Test Strips (no refrigeration required); bar-code protectors; patient and pin-feed labels; storage bags & racks; vapor detectors; a complete line of products for formaldehyde reprocessing; and much more.

### For Product Catalog, Orders or Technical Support:

**Call: 1-800-647-3873, Fax: 1-877-352-5557 or Email: [orders@rpc-rabrenco.com](mailto:orders@rpc-rabrenco.com)**

Part Number	Description	Packaging
MX-30004	4 gallon containers per case	3 liters in each cont.
K100-0100	Micro-X® Residual Test Strips	100 per bottle
K100-0105	Micro-X® Potency Test Strips	100 per bottle